

www.pathoview.com




Pathoview

Personal Storage of Slides



Pathoview

Personal storage for microscopy digital slides



www.pathoview.com

Upload, store and share your digital slides

Communicate and consult with colleagues all over the world



Uploading

Upload your slides in any way convenient for you: via a web browser, from the Vision and Vision Slide Viewer applications, or send them using your Vision Suite organization.



Storage

Save digital slides in Pathoview, and they will always be at your fingertips, wherever you are. Ensure access to digital slides from any device.



Shared access

Information exchange — share and discuss your slides with colleagues. Quick and easy sharing of links to slides of any size.



Digital slides

The service supports the major digital slide formats: .mrxs, .svs, .tiff, .ndpi, .dcm.

More ways to work with digital slides



Glass slide viewing mode






Table view mode




Uploading slides via browser



Creating a link to send to your colleagues for discussion




Digital slide viewing with low resolution




Digital slide viewing with high resolution

Vision and Pathoview integration



In Vision select a sample and click the “Submit” button on the toolbar




In Pathoview click the “Submit” button, copy the link and send it to your colleagues for consultation



Your colleague views the digital slide on any device

Vision Ecosystem Integration



Pre-classification has to be validated by the specialist reviewing the galleries of cells.
 Please check the intended use and the registration in your local regulation.
 Depending on the regulations of the respective region, some products may be used for research purpose only.
 We reserve the right to change specification without notice.