


Digital Microscopy for Diagnostics of Sexually Transmitted Diseases



Application Module Vision Cyto STD (Sexually Transmitted Diseases)




Automatic scanning of slides

Digital slide preparation


Gallery presentation of:

- Trichomonas vaginalis
- Neisseria gonorrhoeae
- Candida spp.
- Clue cells



Trichomonas vaginalis

Candida Spp.



Neisseria gonorrhoeae

Clue cells

Testing for sexually transmitted infections.

Administrative Modules Vision Manager

		Code	Modified by
7 S			WESTMEDICA\o.ku
۰ 🗆			


Analysis procedure automation, data processing rules.

Consulting and Education Modules Vision Suite Vision Expertise

	•• <				vision	-suite.com		Ċ	+
0	WEST Samples (138)						West M	Aedica 🚍 👯 📵
C	T								i
	Slide View		0 U	IS Sample Date	ID	Sample	Size	Note	
	•	MASSITIAN	0	12/05/2018, 5:49 PM	2850	Pap	1.1 GB		
	0	25004500000	0	12/05/2018, 5:46 PM	2849	Pap	437.4 MB		
		State of the second	0	12/05/2018, 5:44 PM	2848	Pap	383.3 MB		
	0	85 9191	0	12/05/2018, 5:39 PM	2847	Рар	403.3 MB		
		TTI NAMAN		12/05/2018, 7:33 PM	4914	DS	1.9 GB		
	102	1-14563-1 1000 - 100 - 100		12/05/2018, 7:33 PM	4558	Histo	905.7 MB		
	0.0			12/05/2018, 7:32 PM	4557	Histo	1.5 GB		

Cloud/server for telemedicine and remote consultations with colleagues.

Vision Remote



Remote workstation.

VISION Vision-sxpertise.com VISION WBC results Q199b 796b Tomas techands Mais tempon Meterson Mais tempon Meterson Mais tempon

Online testing and quality control.

Artificial intelligence



Artificial intelligence (Al) cumulates algorithms and technologies that allow computers to learn and solve intellectual tasks provided by humans.

Al speeds up processing and interpretation of data, and allows to efficiently perform the most comprehensive tasks, including medical image analysis.

Clinical applications



The latest developments of artificial intelligence provide a solution for the tasks connected with automation in digital microscopy.

Our technologies speed up the diagnostic process, reduce analysis time and lower subjectivity of the results received.

They improve the efficiency of laboratory routine operation, bringing microscopy analyses in line with state-of-the-art standards.

Specifications



Vision Assist Cell Imaging Analyzer



Vision Pro Cell Imaging Analyzer



Vision Ultimate Cell Imaging Analyzer

Application module: Vision Cyto STD	Application module: Vision Cyto STD	Application module: Vision Cyto STD
Working modes: queue (only 4 slides version)	Working modes: queue, random access	Working modes: queue, continuous and random access, STAT testing, 24/7
Automatic scanning	Automatic scanning	Automatic scanning
1 or 4 slides	8 slides	Up to 200 slides
Manual slide handling	2 cassettes with slides	Automatic slide handling
Microscope for scanning	Microscope for scanning	Microscope for scanning
Personal computer	Personal computer	Personal computer
Monitor	Monitor	Monitor
—	—	Touch-screen monitor for control
-	Built-in barcode reader (optional)	Built-in barcode reader
_	Automatic oil dispenser (optional)	Automatic oil dispenser
Optical system: 10x, 50x Oil, 100x Oil	Optical system: 10x, 60x Oil, 100x Oil	Optical system: 10x, 50x Oil, 100x Oil
Bright field	Bright field	Bright field
Köhler illumination, LED	Köhler illumination, LED	Köhler illumination, LED
Bidirectional LIS, LIS2-A2 (ASTM), HL7, Ethernet	Bidirectional LIS, LIS2-A2 (ASTM), HL7, Ethernet	Bidirectional LIS, LIS2-A2 (ASTM), HL7, Ethernet

Art. N.: 71150.00.019 (1 slide)	Art. N.: 72852.00.019 (8 slides)	Art. N.: 73011.00.019 (200 slides)
Art. N.: 71450.00.019 (4 slides)		

We reserve the right to change specification without notice.



West Medica Produktions- und Handels-GmbH Brown-Boveri-Straße 6, B17-1 2351 Wiener Neudorf, Austria tel.: +43 (0) 2236 892465, fax: +43 (0) 2236 892464 vienna@westmedica.com, www.wm-vision.com

Official Distributor