


www.pathoview.com


Pathoview

Upload, Store and Share your Digital Slides

Communicate and Consult with Colleagues all over the World



www.pathoview.com



Your Personal Storage for Microscope Digital Slides



Uploading

Upload your slides in any way convenient for you: via a web browser, from the Vision and Vision Slide Viewer applications, or send them using your Vision Suite organization.



Storage

Save digital slides in Pathoview, and they will always be at your fingertips, wherever you are. Ensure access to digital slides from any device.



Shared access

Information exchange — share and discuss your slides with colleagues. Quick and easy sharing of links to slides of any size.



Digital slides

The service supports the major digital slide formats: .mrxs, .svs, .tiff, .ndpi, .dcm.

Support for Slide Scanners from Various Manufacturers



West Medica (.tiff)



3DHISTECH (.mrxs)



Leica (.svs)




Hamamatsu (.ndpi)



Philips (.isyntax)



DICOM (.dcm)



Glass slide view mode







Table view mode




Uploading slides via browser



Creating a link to send to your colleagues for discussion





Digital slide viewing with low magnification



Digital slide viewing with high magnification


More Ways to Work with Digital Slides






Histopathology

Cytology






Hematology

Microbiology

Vision Products Integration



We reserve the right to change specification without notice.

